

CERTIFICATE OF ANALYSIS FOR  
**GOLD ORE REFERENCE MATERIAL**  
**OREAS 17Pb**

SUMMARY STATISTICS

Constituent	Recommended Value	95% Confidence Interval		Tolerance limits 1- $\alpha$ =0.99, $\rho$ =0.95	
		Low	High	Low	High
Gold, Au (ppm)	2.56	2.50	2.62	2.54	2.58

*Prepared by:*  
*Ore Research & Exploration Pty Ltd*  
*August, 2003*

## INTRODUCTION

OREAS certified reference materials (CRMs) are intended to provide a low cost method of evaluating and improving the quality of precious and base metal analysis of geological samples. To the analyst they provide an effective means of calibrating analytical equipment, assessing new techniques and routinely monitoring in-house procedures. To the geologist they provide a means of implementing quality control in analytical data sets generated in exploration, from the grass roots level through to prospect evaluation, and in grade control at mining operations.

As a rule only source materials exhibiting a high level of homogeneity of the element(s) of interest are used in the preparation of these materials. This has enabled Ore Research & Exploration to produce a range of gold ore CRMs exhibiting homogeneity that matches or exceeds that of currently available international reference materials. In certain instances CRMs produced from a single source are sufficiently homogeneous to produce a relatively coarse-grained form designed to simulate drill chip samples. These have a grain size of minus 3mm and are designated with a "C" suffix to the CRM identification number. These standards are packaged in 1kg units following homogenisation and are intended for submission to analytical laboratories in subsample sizes of as little as 250g. They offer the added advantages of providing a check on both sample preparation and analytical procedures while acting as a blind standard to the assay laboratory. The more conventional pulped standards have a grain size of minus 20 to minus 75 microns and a higher degree of homogeneity. These standards are distinguished by a "P" suffix to the standard identification number. In line with ISO recommendations successive batch numbers are now designated by the lower case suffixes "a", "b", "c", "d", etc.

## SOURCE MATERIALS

Reference material OREAS 17Pb was prepared from a blend of barren alkali olivine basalt from Epping, Victoria, Australia and gold-bearing Magdala ore from the Stawell Gold Mine, west-central Victoria, Australia. The Magdala lode is intimately associated with an intensely deformed package of volcanogenic sedimentary rocks. Mineralisation in the ore consists of a quartz-sericite-carbonate schist assemblage containing the sulphides pyrite and arsenopyrite. The major constituents of the alkali olivine basalt are feldspar, augite, olivine and titanomagnetite.

The approximate major and trace element composition of this sulphide-bearing gold ore standard OREAS 17Pb is given in Table 1. The constituents SiO<sub>2</sub> to Total and Ba, Ni, V and Zr are the means of duplicate XRF analyses determined using a borate fusion method, C and S are the means of duplicate Leco analyses, while the remaining constituents, As to Yb, are means of twenty replicate analyses determined via INAA at Becquerel Laboratories.

Gold homogeneity has been evaluated and confirmed by INAA on twenty-four 0.5 gram sample portions and by a nested ANOVA program using conventional fire assay. The tolerance interval is determined from the INAA data while the recommended value and confidence interval are based on a round robin program incorporating a total of 76 analyses at 17 laboratories.

## COMMUNITION AND HOMOGENISATION PROCEDURES

The gold-bearing basaltic material comprising OREAS 17Pb was prepared in the following manner:

- a) *jaw crushing to minus 7mm*
- b) *drying to constant mass at 105°C*
- c) *milling of the barren basalt to 98% minus 75 micron*
- d) *milling of the gold ore to 100% minus 20 micron*
- e) *blending in appropriate proportions to achieve the desired grade*
- f) *bagging into 25kg sublots*

Throughout the bagging stage twenty-four 1kg test units were taken at regular intervals, sealed in laminated plastic bags and set aside for the analytical program.

Table 1. Indicative major and trace element composition of gold ore reference material OREAS 17Pb; SiO<sub>2</sub> to Total as weight percent; rest in parts per million; SiO<sub>2</sub> to Zr by fusion XRF except C and S by Leco furnace; As to Yb by INAA.

Constituent	Concentration (XRF)	Constituent	Concentration (INAA)
SiO <sub>2</sub>	53.75	As	3933
TiO <sub>2</sub>	1.46	Ce	35
Al <sub>2</sub> O <sub>3</sub>	12.4	Co	40
Fe <sub>2</sub> O <sub>3</sub>	12.8	Cr	224
MnO	0.21	Cs	<2
MgO	6.16	Eu	1.4
CaO	7.34	Hf	2
Na <sub>2</sub> O	2.51	La	16
K <sub>2</sub> O	0.77	Lu	<0.5
P <sub>2</sub> O <sub>5</sub>	0.33	Rb	<40
LOI	1.66	Sb	2.9
<b>Total</b>	<b>99.39</b>	Sc	17
Ba	280	Sm	4.5
C (Leco)	0.255	Sr	0.04
Ni	215	Th	<5
S (Leco)	1.34	W	<4
V	265	Yb	1.5
Zr	142.5		

## ANALYSIS OF OREAS 17Pb

At the time of writing seventeen laboratories had reported their results and these respondents are listed in the section headed Participating Laboratories.

To maintain anonymity laboratories have been randomly designated the letter codes A through Q. With the exception of Laboratory Q, each laboratory received four 110g samples with instructions to carry out one 30 to 50g fire assay determination for gold on each sample. Apart from Lab C (ICPOES finish) and Lab G (gravimetric finish) the remaining laboratories employed a flame AAS finish.

For each laboratory two 110g subsamples were scoop-split from each of two separate 1kg test units taken during the bagging stage. This two-stage nested design for the

interlaboratory programme was amenable to analysis of variance (ANOVA) treatment and enabled a comparative assessment of within- and between-unit homogeneity.

For the determination of a statistical tolerance interval, a 10g scoop split was taken from each of the twenty-four random test units and submitted to Laboratory Q for gold assay via instrumental neutron activation analysis on a reduced analytical subsample weight of 0.5g.

Individual assay results for the fire assay and INAA methods are presented in Tables 2 and 3 together with the mean, median and standard deviation (absolute and relative) for each data set. Interlaboratory agreement of the means is good for all but two data sets, Labs A and O, where bias is unacceptably high (+16.6%) and low (-10.4%), respectively.

## STATISTICAL EVALUATION OF ANALYTICAL DATA FOR OREAS 17Pb

### Recommended Value and Confidence Limits

The recommended value was determined from the mean of means of accepted replicate values of accepted laboratory data sets A to K and N to P according to the formulae

$$\bar{x}_i = \frac{1}{n_i} \sum_{j=1}^{n_i} x_{ij}$$

$$\bar{x} = \frac{1}{p} \sum_{i=1}^p \bar{x}_i$$

where

$x_{ij}$  is the  $j$ th result reported by laboratory  $i$ ;

$p$  is the number of participating laboratories;

$n_i$  is the number of results reported by laboratory  $i$ ;

$\bar{x}_i$  is the mean for laboratory  $i$ ;

$\bar{x}$  is the mean of means.

The confidence limits were obtained by calculation of the variance of the consensus value (mean of means) and reference to Student's- $t$  distribution with degrees of freedom ( $p-1$ )

$$\hat{V}(\bar{x}) = \frac{1}{p(p-1)} \sum_{i=1}^p (\bar{x}_i - \bar{x})^2$$

$$\text{Confidence limits} = \bar{x} \pm t_{1-x/2}(p-1) \left( \hat{V}(\bar{x}) \right)^{1/2}$$

where  $t_{1-x/2}(p-1)$  is the  $1-x/2$  fractile of the  $t$ -distribution with  $(p-1)$  degrees of freedom.

The distribution of the values are assumed to be symmetrical about the mean in the calculation of the confidence limits.

Table 2. Analytical results for gold in OREAS 17Pb (INAA - instrumental neutron activation analysis; FA\*AAS - fire assay / atomic absorption spectrometry; FA\*ICP - fire assay/inductively coupled plasma optical emission spectrometry; FA\*GRAV - fire assay/gravimetric; Std. Dev. and Rel.Std.Dev. are one sigma values; PDM<sup>3</sup>- percent deviation of lab mean from corrected mean of means; outliers in bold; values in ppm).

Sample No.	Lab A INAA (30g)	Lab B FA*AAS (50g)	Lab C FA*ICP (40g)	Lab D FA*AAS (50g)	Lab E FA*AAS (50g)	Lab F FA*AAS (50g)	Lab G FA*GRAV (50g)	Lab H FA*AAS (50g)	Lab I FA*AAS (50g)	Lab J FA*AAS (50g)	Lab K FA*AAS (50g)
1	2.94	2.59	2.43	2.59	2.54	2.47	2.68	2.58	2.58	2.63	2.76
2	3.02	2.56	2.62	2.57	2.60	2.58	2.69	2.56	2.47	2.72	2.65
3	2.98	2.66	<b>1.96</b>	2.58	2.62	2.54	2.67	2.58	2.54	2.70	
4	2.99	2.60	2.46	2.58	2.64	2.49	2.68	2.57	2.52	2.69	
Mean:	<b>2.98</b>	2.60	2.37	2.58	2.60	2.52	2.68	2.57	2.53	2.69	2.70
Median:	2.99	2.60	2.45	2.58	2.61	2.52	2.68	2.58	2.53	2.70	2.70
STDEV:	0.03	0.04	0.28	0.01	0.04	0.05	0.01	0.01	0.05	0.04	0.08
RSD:	1.11%	1.61%	12.00%	0.32%	1.66%	1.97%	0.31%	0.37%	1.81%	1.44%	3.09%

Table 2. continued

Sample No.	Lab L FA*AAS (50g)	Lab M FA*AAS (50g)	Lab N FA*AAS (50g)	Lab O FA*AAS (50g)	Lab P FA*AAS (30g)
1	2.60	2.27	2.52	2.37	2.38
2	2.65	2.36	2.68	2.44	2.21
Mean:	2.63	2.32	2.60	2.41	<b>2.29</b>
Median:	2.63	2.32	2.60	2.41	2.29
STDEV:	0.04	0.06	0.12	0.05	0.12
RSD:	1.35%	2.75%	4.52%	5.21%	2.06%

Table 3. Analytical results for gold (ppm) in OREAS 17Pb by instrumental neutron activation analysis on 0.5g analytical subsample weights (abbreviations as for Table 2).

Sample No.	Lab Q
1	2.55
2	2.46
3	2.52
4	2.47
5	2.43
6	2.47
7	2.32
8	2.54
9	2.40
10	2.47
11	2.51
12	2.44
13	2.52
14	2.44
15	2.34
16	2.40
17	2.44
18	2.48
19	2.54
20	2.50
21	2.41
22	2.33
23	2.40
24	<b>2.27</b>
Mean:	2.44
Median:	2.45
STDEV:	0.08
RSD:	3.10%

The test for rejection of individual outliers was based on the test criterion, T, and reference to tables of critical values of T at the 1% level of significance (ASTM E 178-94) as follows:

$$T_{ij} = \left| (x_{ij} - \bar{x}_i) \right| / s_i$$

where

$T_{ij}$  is the test criterion for the  $j$ th result of laboratory  $i$ ;  
 $s_i$  is the standard deviation of laboratory  $i$ .

The same principles were applied in testing for outlying laboratory means. Individual and mean outliers are shown in bold type in Tables 2 and 3 and have been omitted in the determination of recommended values.

Table 4. Recommended value and 95% confidence interval

Constituent	Recommended value	95% Confidence interval	
		Low	High
Gold, Au (ppm)	2.56	2.50	2.62

### Statement of Homogeneity

The variability of replicate assays from each laboratory is a result of both measurement and subsampling errors. In the determination of a statistical tolerance interval it is therefore necessary to eliminate, or at least substantially minimise, those errors attributable to measurement. One way of achieving this is by substantially reducing the analytical subsample weight to a point where most of the variability in replicate assays is due to inhomogeneity of the reference material and measurement error becomes negligible. This approach was adopted in the INAA data set (Table 3) where a 0.5 g subsample weight was employed. The homogeneity was determined from tables of factors for two-sided tolerance limits for normal distributions (ISO Guide 3207) in which

$$\text{Lower limit is } \bar{x} - k'_2(n, p, 1 - \alpha)s$$

$$\text{Upper limit is } \bar{x} + k'_2(n, p, 1 - \alpha)s$$

where

$n$  is the number of results reported by laboratory  $Q$ ;

$1 - \alpha$  is the confidence level;

$p$  is the proportion of results expected within the tolerance limits;

$k'_2$  is the factor for two-sided tolerance limits ( $m, \sigma$  unknown);

and  $s$  is computed according to the formula

$$s = \left[ \frac{\sum_{j=1}^n (x_j - \bar{x})^2}{n-1} \right]^{1/2}$$

**No individual outliers were removed from the results prior to the calculation of tolerance intervals.**

From the INAA data set an estimated tolerance interval of  $\pm 0.02$  ppm at an analytical subsample weight of 50g was obtained (using the sampling constant relationship of Ingamells and Switzer, (1973) and is considered to reflect the actual homogeneity of the material under test. The meaning of this tolerance interval may be illustrated for gold (refer Table 5), where 99% of the time at least 95% of 50g-sized subsamples will have concentrations lying between 2.54 and 2.58 ppm. Put more precisely, this means that if the same number of subsamples were taken and analysed in the same manner repeatedly, 99% of the tolerance intervals so constructed would cover at least 95% of the total population, and 1% of the tolerance intervals would cover less than 95% of the total population (ISO Guide 35).

Table 5. Recommended value and tolerance interval.

Constituent	Recommended value	Tolerance interval $1-\alpha=0.99, \rho=0.95$	
		Low	High
Gold, Au (ppm)	2.56	2.54	2.58

## PARTICIPATING LABORATORIES

Acme Analytical Laboratories Ltd, Vancouver, Ontario, Canada  
 ALS Chemex Laboratories Pty Ltd, North Vancouver, Ontario, Canada  
 ALS Chemex Laboratories Pty Ltd, Garbutt, Queensland, Australia  
 ALS Chemex Laboratories Pty Ltd, Sparks, Nevada, USA  
 Amdel Laboratories, Wangara, Western Australia, Australia  
 Amdel Laboratories, Thebarton, South Australia, Australia  
 SGS Analabs, Welshpool, Western Australia, Australia  
 Becquerel Laboratories, Lucas Heights, New South Wales, Australia  
 Becquerel Laboratories, Mississauga, Ontario, Canada  
 Genalysis Laboratory Services Pty Ltd, Maddington, WA, Australia  
 Indo Assay Laboratories, Kalimantan Timur, Indonesia  
 Intertek Utama Services, Jakarta Selatan, Indonesia  
 SGS Lakefield Research, Lakefield, Ontario, Canada  
 Lakefield Geosol Ltda, Belo Horizonte - MG, Brazil  
 McPhar Geoservices (Phil) Inc., Makati, Manilla, Philippines  
 OMAC Laboratories, Loughrea, County Galway, Ireland  
 Standard & Reference Laboratories, Malaga, Western Australia, Australia

## **PREPARER AND SUPPLIER OF THE REFERENCE MATERIAL**

The gold ore reference material, OREAS 17Pb has been prepared and certified and is supplied by:

Ore Research & Exploration Pty Ltd  
3 London Drive  
Bayswater VIC 3153  
AUSTRALIA

Telephone	(03) 9762 1808	International	+613-9762 1808
Facsimile	(03) 9762 3808	International	+613-9762 3808
Email	info@ore.com.au	Web	www.ore.com.au

It is available in unit sizes of 60g laminated foil packets.

## **INTENDED USE**

OREAS 17Pb is a reference material intended for the following:

- i) for the calibration of instruments used in the determination of the concentration of gold;
- ii) for the verification of analytical methods for gold;
- iii) for the preparation of secondary reference materials of similar composition;
- iv) as an arbitration sample for commercial transactions.

## **STABILITY AND STORAGE INSTRUCTIONS**

OREAS 17Pb has been prepared from a blend of gold-ore and barren basalt. Being characterised by a low sulphide content the sealed CRM is considered to have long-term stability under normal storage conditions.

## **INSTRUCTIONS FOR THE CORRECT USE OF THE REFERENCE MATERIAL**

The recommended value for OREAS 17Pb refers to the concentration level of gold in its packaged form at hygroscopic equilibrium. Drying in air to constant mass at 105<sup>0</sup>C has established a hygroscopic moisture content of 1.20%. If the reference material is dried by the user prior to analysis, the recommended value should be corrected to the moisture-free basis.

## **LEGAL NOTICE**

Ore Research & Exploration Pty Ltd has prepared and statistically evaluated the property values of this reference material to the best of its ability. The Purchaser by receipt hereof releases and indemnifies Ore Research & Exploration Pty Ltd from and against all liability and costs arising from the use of this material and information.

**CERTIFYING OFFICER:** Dr Paul Hamlyn



## **ACKNOWLEDGMENTS**

The cooperation of all participating laboratories is gratefully acknowledged.

## **REFERENCES**

ASTM E178-94 (1994), Standard practice for dealing with outlying observations.

Ingamells, C. O. and Switzer, P. (1973), *Talanta* 20, 547-568.

ISO Guide 35 (1985), Certification of reference materials - General and statistical principals.

ISO Guide 3207 (1975), Statistical interpretation of data - Determination of a statistical tolerance interval.